

**ABDOS** 

 **ABDOS**   
LIFE SCIENCE



# ABDOS Pipette Tips

Creating For Life Sciences

# THE MOST COMPATIBLE PIPETTING EXPERIENCE WITH ABDOS TIPS



- Abdos pipette tips are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & confirming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- The use of the latest injection moulding processes produces high precision, thin walled tips with high clarity for better visibility of samples and ensuring perfect seal while reducing the loading and ejection force.
- The universal tip collar fits all popular pipettes.
- Diamond polished moulds for an extremely hydrophobic & smooth surface to ensure lowest retention in its class.
- Perfect air lock grip on pipette by inner radial rings for accurate pipetting.
- The use of multiple cavity tools guarantee high conformity within one batch. The automated moulding and assembling procedures carried out under validated clean-room conditions to eliminate contamination.



- Each wafer of the rack is printed with the effective vol range and lot no. with product code. It makes easy to keep your work organized and ensure supplies are traceable.
- 4 stage quality control procedures ensures 100% quality & performance.
- Sterile tips has been validated according to ENISO 11137-2: 2012 for SAL  $10^{-6}$ .



FREE



FREE



FREE



FREE



FREE



121 °C

