

**ABDOS** 

 **ABDOS**   
LIFE SCIENCE



# ABDOS CENTRIFUGE TUBES

# Experience the Consistency & Purity with Abdos Centrifuge Tubes

Confirms US FDA 21 CFR, **free from natural rubber, heavy metal & inhibitory plasticizers.**


**Even wall thickness** to withstand maximum a force, 20,000 RCF.

Sterile tube have been validated according to **ENISO 11137-2: SAL 10-6.**

**Precise lid sealing & perfect grip on screw** for secure control on the cap in wet conditions.

**Highly Transparent -** 5% less haze percentage than a competitor.

**Suitable for temperature of -80° to +121° C**



**Larger white frosted portion** for easy writing.

**Robust polymer selection -** tubes do not leech & migrate unwanted substances into your sample.

Unique printing technology - chemical resistance, scratch proof for **accurate & legible volume measurement.**

**Hydrodynamic Conical bottom** for perfect peleting of sediments.

## The Caps

**Flat thread confirm vacuum tight closure -** leak proof and prevent from override of cap.



**Unique plug seal -** ensuring perfect sealing & make it leak-proof.