


Filter Pipette Tips

Creating For Life Sciences

Abdos Super Filter Prevents Both Liquid & Aerosol Passing Through The Filter.

Competitor Filter pipette Tips

Abdos Super Filter Pipette Tips




Abdos Super Filter Technology provides aerosol barrier & aqueous liquid* by-pass prevention in over-pipetting thus safeguarding precious samples against contamination.

Ensure Smooth & Hassle Free Pipetting with Increased Air Gap to Prevent Cross Contamination.

Abdos 200 µl Filter Tips

INCREASED AIR GAP

Competitor 200 µl Filter Tips



Experience The Finest Super Filter Technology With Abdos Filter Tips



- Universal fit, leak proof fitment.
- Precise safety gap between filter and liquid, filter and pipette shaft to avoid contamination & worry free pipetting.
- Abdos Filter is 100% pure porous polyethylene (PE) helps us to recover samples from wetted filter and does not trap or block the liquid into the filter.
- 5 Stage quality control procedures ensures 100% quality & performance.
- Inert hydrophobic filter does not trap or block the liquid, thus save your precious sample loss in pipetting error.
- Sterile tips has been validated according to ENISO 11137-2; 2012 for SAL 10^{-6} .
- Abdos filter pipette tips are manufactured from ultra pure USP class VI certified virgin polypropylene to ensure the highest quality & confirming to US FDA 21 CFR. Free from heavy metals, natural rubber & inhibitory plasticizers.
- Proprietary polyethylene super filter provides aerosol barrier and aqueous liquid by-pass prevention in over-pipetting hereby preventing cross contamination.
- Best suited for biological application.